


VACUMED® SODIO CITRATO 3,2 %
 VACUMED® SODIUM CITRATE 3,2 %

VACUMED® SODIO CITRATO 3,8 %
 VACUMED® SODIUM CITRATE 3,8 %



42410

42411

42412

42480

42481

42482

Vacumed®
 13x75 mm con Sodio
 Citrato 3,2 % x 1,8 ml
 di sangue per
 coagulazione, tappo
 azzurro, sterile

Vacumed®
 13x75 mm con Sodio
 Citrato 3,2 % x 2,7 ml
 di sangue per
 coagulazione, tappo
 azzurro, sterile

Vacumed®
 13x75 mm con Sodio
 Citrato 3,2 % x 3,6 ml
 di sangue per
 coagulazione, tappo
 azzurro, sterile

Vacumed®
 13x75 mm con Sodio
 Citrato 3,8 % x 1,8 ml
 di sangue per
 coagulazione, tappo
 azzurro, sterile

Vacumed®
 13x75 mm con Sodio
 Citrato 3,8 % x 2,7 ml
 di sangue per
 coagulazione, tappo
 azzurro, sterile

Vacumed®
 13x75 mm con Sodio
 Citrato 3,8 % x 3,6 ml
 di sangue per
 coagulazione, tappo
 azzurro, sterile

Vacumed® 13x75 mm
 with Sodium Citrate
 3,2 % x 1,8 ml
 of blood for
 coagulation, light
 blue cap, sterile

Vacumed® 13x75 mm
 with Sodium Citrate
 3,2 % x 2,7 ml
 of blood for
 coagulation, light
 blue cap, sterile

Vacumed® 13x75 mm
 with Sodium Citrate
 3,2 % x 3,6 ml
 of blood for
 coagulation, light
 blue cap, sterile

Vacumed® 13x75 mm
 with Sodium Citrate
 3,8 % x 1,8 ml
 of blood for
 coagulation, light
 blue cap, sterile

Vacumed® 13x75 mm
 with Sodium Citrate
 3,8 % x 2,7 ml
 of blood for
 coagulation, light
 blue cap, sterile

Vacumed® 13x75 mm
 with Sodium Citrate
 3,8 % x 3,6 ml
 of blood for
 coagulation, light
 blue cap, sterile