



TEST TUBE CYLINDRICAL SORWALL TYPE

tubo cilíndrico tipo Sorwall / tube cylindrique type Sorwall / provetta cilindrica tipo Sorwall



EXCELLENT RESISTANCE

excelente resistencia / excellente résistance / eccellente resistenza



POLYPROPYLENE

polipropileno / polypropylène / polipropilene

	CODE	CAPACITY	TYPE	PACKAGING
	21014	5 ml	12x75 mm	IN BAG OF 250 PCS en bolsa de 250 piezas / en sachet de 250 pièces / in sacchetto da 250 pz
	21014B	5 ml	12x75 mm	IN BULK a granel / en vrac / sfusa