



FAECES CONTAINERS 30 ML

Frascos 30 ml para heces / Pots 30 ml pour selles / Contenitori feci 30 ml



FAECES CONTAINER WITH FROSTED LABEL

frasco para heces con banda mate / pot pour selles avec surface d'écriture / contenitore per feci con superficie di scrittura

WITH SCREW CAP APART

con tapón no roscado / avec bouchon à vis no jointé / con tappo a vite a parte




POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène / polipropilene

	CODE	CAPACITY	STERILITY	PACKAGING
	25130	30 ml (white)	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa
	25131	30 ml (red)	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa