



# STERILE TEST TUBES

Tubos de ensayo estériles / Tubes stériles / Provette sterili



## TEST TUBE CYLINDRICAL

tubo cilíndrico / tube cylindrique /  
provetta cilindrica

## MAXIMUM TRANSPARENCY

máxima transparencia / transparence maximale /  
massima trasparenza

## EXCELLENT RESISTANCE

excelente resistencia / excellente résistance /  
eccellente resistenza

## STERILE

estéril / stérile / sterile

## RED

rojo / rouge / rosso





## POLYETHYLENE

polietileno / polyéthylène / polietilene



## POLYSTYRENE

poliestireno / polystyrène / polistirolo

|                                                                                     | CODE  | CAPACITY | TYPE                       | PACKAGING                                                                                    |
|-------------------------------------------------------------------------------------|-------|----------|----------------------------|----------------------------------------------------------------------------------------------|
|  | 25091 | 12 ml    | 16x100 mm<br>(cylindrical) | INDIVIDUALLY WRAPPED<br>embalado individualmente / sachet<br>unitaire / in sacchetto singolo |
|  | 25092 | 12 ml    | 16x100 mm<br>(conical)     | INDIVIDUALLY WRAPPED<br>embalado individualmente / sachet<br>unitaire / in sacchetto singolo |