


PRE-ASSEMBLED TEST TUBE RACKS

Gradillas para tubos premontadas / Racks pour tubes prèmontès /
Portaprovette premontati




RACK FOR TEST TUBES

gradilla para tubos / portoir pour tubes /
portaprovette



AUTOCLAVABLE

autoclavable / autoclavable / autoclavabile



POLYPROPYLENE

polipropileno / polypropylène /
polipropilene

	CODE	COMPATIBILITY	CAPACITY
	27152	test tubes Ø 12 and 13 mm	90 places
	27156	test tubes Ø 16 mm	60 places