



**TEST TUBE WITH LITHIUM HEPARIN**

tubo con Heparina de Litio / tube avec Héparine de Lithium / provetta con Litio Eparina

**PERFECT TRACEABILITY**

trazabilidad perfecta / traçabilité parfaite / tracciabilità perfetta

**UNBREAKABLE**

irrompibile / incassable / infrangibile

**DARK GREEN**

verde oscuro / vert foncé / verde scuro




**POLYETHYLENE**

polietileno / polyéthylène / polietilene



**POLYPROPYLENE**

polipropileno / polypropylène / polipropilene

	CODE	CAPACITY	TYPE	PACKAGING
	32306	2,5 ml	12x56 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32312	2,5 ml	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32314	3 ml	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32316	4 ml	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32304	5 ml	12x86 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi