



TEST TUBE WITH 0,25 ML OF SODIUM CITRATE (3,8%) X 1 ML OF BLOOD (ESR) AND 2,25 ML OF BLOOD (COAGULATION)

tubo con 0,25 ml de Citrato de Sodio (3,8%) x 1 ml de sangre (VSG) y 2,25 ml de sangre (Coagulación) / tube avec 0,25 ml de Citrate de Sodium (3,8 %) x 1 ml de sang (VS) y 2,25 ml de sang (Coagulation) / provetta con 0,25 ml di Sodio Citrato (3,8%) x 1 ml di sangue (VES) e 2,25 ml di sangue (Coagulazione)

PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite / tracciabilità perfetta

UNBREAKABLE

irrompible / incassable / infrangibile

PINK

rosa / rose / rosa




POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène / polipropilene

	CODE	CAPACITY	TYPE	PACKAGING
	22213	1 ml x ESR 2,25 ml x COAGULATION	12x56 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22218	1 ml x ESR 2,25 ml x COAGULATION	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22210	1 ml x ESR 2,25 ml x COAGULATION	12x86 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi