




# SODIUM CITRATE TUBES FOR ESR

Tubos con Citrato de Sodio para VSG / Tubes avec Citrate de Sodium pour VS /  
Provette con Sodio Citrato per VES



## TEST TUBE WITH 0,4 ML OF SODIUM CITRATE (3,8%) X 1,6 ML OF BLOOD (ESR)

tubo con 0,4 ml de Citrato de Sodio (3,8%) x 1,6 ml de sangre (VSG) / tube avec 0,4 ml de Citrate de Sodium (3,8 %) x 1,6 ml de sang (VS) / provetta con 0,4 ml di Sodio Citrato (3,8%) x 1,6 ml di sangue (VES)

### PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite /  
tracciabilità perfetta

### UNBREAKABLE

irrompibile / incassable / infrangibile



### POLYETHYLENE

polietileno / polyéthylène / polietilene



### POLYPROPYLENE

polipropileno / polypropylène /  
polipropilene

	CODE	CAPACITY	TYPE	PACKAGING
	22212	1,6 ml x ESR	12x56 mm (pink)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22257	1,6 ml x ESR	13x75 mm (pink)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22209	1,6 ml x ESR	12x86 mm (pink)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32267	1,6 ml x ESR	13x75 mm (black)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi