



**TEST TUBE WITH 0,2 ML OF SODIUM CITRATE (3,2%) X 1,8 ML OF BLOOD (COAGULATION), PIERCEABLE CAP**

tubo con 0,2 ml de Citrato de Sodio (3,2%) x 1,8 ml de sangre (Coagulación), tapón perforable / tube avec 0,2 ml de Citrate de Sodium (3,2%) x 1,8 ml de sang (Coagulation), bouchon perforable / provetta con 0,2 ml di Sodio Citrato (3,2%) x 1,8 ml di sangue (Coagulazione), tappo perforabile

**PERFECT TRACEABILITY**

trazabilidad perfecta / traçabilité parfaite / tracciabilità perfetta

**UNBREAKABLE**

irrompible / incassable / infrangibile

**LIGHT BLUE**

azul claro / bleu clair / azzurro



**RUBBER**

goma / caoutchouc / gomma



**POLYPROPYLENE**

polipropileno / polypropylène / polipropilene

CODE	CAPACITY	TYPE	PACKAGING
32145	0,2 ml of Sodium x 1,8 ml of blood	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi