



TEST TUBE WITH CLOT ACTIVATOR SILICONE COATED INTERIOR

tubo con activador de la coagulación interior recubierto de silicona / tube avec activateur de la coagulation silicone à l'intérieur / provetta con attivatore della coagulazione siliconata internamente

PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite / tracciabilità perfetta

UNBREAKABLE

irrompibile / incassable / infrangibile

RED

rojo / rouge / rosso





POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène / polipropilene

	CODE	TYPE	PACKAGING
	22942	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22922	16x100 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi