



**PETRI DISH WITHOUT VENTS**

placa de petri sin vientos / boîte de petri sans ergots / piastra di petri senza ventilazione

**MAXIMUM TRANSPARENCY**

máxima transparencia / transparence maximale / massima trasparenza

**PERFECT PLANARITY**

planaridad perfecta / parfaite planéité / planarità perfetta



**POLYSTYRENE**

poliestireno / polystyrène / polistirolo

	CODE	TYPE	STERILITY	PACKAGING
	29142	ø 90 mm, H 14,2 mm	PRODUCED IN CLEANROOM ISO 6 producido en cleanroom ISO 6 / produit en cleanroom ISO 6 / prodotta in camera bianca ISO 6	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi
	29182	ø 90 mm, H 14,2 mm	STERILE estéril / stérile / sterile	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi
	29152	ø 90 mm, H 16,2 mm	PRODUCED IN CLEANROOM ISO 6 producido en cleanroom ISO 6 / produit en cleanroom ISO 6 / prodotta in camera bianca ISO 6	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi
	29162	ø 90 mm, H 16,2 mm	STERILE estéril / stérile / sterile	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi