




PETRI DISHES

Placas de petri / Boîtes de petri / Piastre di petri



PETRI DISH SINGLE VENT

placa de petri con un viento / boîte de petri avec un ergoit / piastra di petri singola ventilazione

MAXIMUM TRANSPARENCY

máxima transparencia / transparence maximale / massima trasparenza

PERFECT PLANARITY

planaridad perfecta / parfaite planéité / planarità perfetta



POLYSTYRENE

poliestireno / polystyrène / polistirolo

	CODE	TYPE	STERILITY	PACKAGING
	29242	ø 90 mm, H 14,2 mm	PRODUCED IN CLEANROOM ISO 6 producido en cleanroom ISO 6 / produit en cleanroom ISO 6 / prodotta in camera bianca ISO 6	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi
	29282	ø 90 mm, H 14,2 mm	STERILE estéril / stérile / sterile	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi
	29252	ø 90 mm, H 16,2 mm	PRODUCED IN CLEANROOM ISO 6 producido en cleanroom ISO 6 / produit en cleanroom ISO 6 / prodotta in camera bianca ISO 6	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi
	29262	ø 90 mm, H 16,2 mm	STERILE estéril / stérile / sterile	BAG OF 20 PIECES bolsa de 20 piezas / en sachet de 20 pièces / in sacchetto da 20 pezzi