

Vacumed® Granules Separators Tubes

FL Medical works to satisfy customers' needs. The company is continuously updating to improve the products and realize new projects.
At the basis of this item there are FL Medical values:

Innovation
Quality
Efficiency

The main characteristics of the Vacumed® Granules Separators Tubes are:

FASTER COAGULATION

The clot activator covers a larger surface of test tube ensuring a faster coagulation

CENTRIFUGE SPEED

The low specific weight of the granules permits less centrifugation force avoiding breakage of blood cells.

PERFECT SERUM SEPARATION

Granules interpose between whole blood and serum obtained after centrifugation



Code	Tube Size	Aspiration Volume	Packaging
42681	13 x 75 mm	3,5 ml	100 pcs per rack 1.200 pcs per box
43686	13 x 100 mm	5 ml	100 pcs per rack 1.200 pcs per box
44688	16 x 100 mm	8 ml	100 pcs per rack 1.200 pcs per box